

LED Technologies, LLC - Clinical Study Information

In a clinical study to determine the efficacy of the DPL® Therapy System, 13 subjects received DPL® Light Therapy treatments over 10 weeks. 100% of the subjects reported visible improvement in fine lines, wrinkles and skin tone.

Reduction of fine lines: 75 % of the subjects rated their improvement between 3 and 5 (on a scale of 5) over the ten week treatment cycle.

Reduction of crow's-feet: 66.6% of the subjects rated their improvement between 3 and 5 (on a scale of 5) over the ten week treatment cycle.

Improvement in skin tone: 91.6% of the subjects rated their improvement between 3 and 5 (on a scale of 5) over the ten week treatment cycle.

Overall skin condition: 91.6% of the subjects rated their improvement between 3 and 5 (on a scale of 5) over the ten week treatment cycle.

Visual Results after 10 Weeks





